INSTRUCTIONS FOR USE

1

Actifuse Microgranules contains phase-pure, porous silicate substituted calcium phosphate granules for use as a synthetic bone graft substitute.

Sterile - Single use - Store in a cool, dry place.

Description:

- 0.8 wt% silicate substituted calcium phosphate
- >99% phase pure silicate substituted calcium phosphate
- Total porosity 80% \pm 2.5%

Indications: as per morsellised grafting.

Directions for use:

- Use in place of or added to morsellised autograft bone.
- Use in place of or added to morsellised allograft bone.
- Actifuse Microgranules may be used alone or mixed with autograft or allograft bone for spinal fusion.
- Actifuse Microgranules may be used alone to fill small unloaded bone defects, e.g. small tumours.
- Actifuse Microgranules used alone may be mixed with either host blood or bone marrow.
- Mixing containers and stirrers must be sterilised before use.
- Actifuse Microgranules implanted alone or when mixed is used in the same manner as morsellised allograft.
- Actifuse Microgranules may be used to fill an inter-body fusion cage.

Risks and warnings:

- The risks with Actifuse Microgranules (e.g. infection, graft loosening, non-union, delayed union, repeat surgery) are the same as those for allograft bone. However, with Actifuse Microgranules there is no risk of disease transmission.
- Actifuse Microgranules should not be used where it could be subject to tension, torsion, compression, shear or bending. A conventional implant (e.g. screw, rod) can protect the graft from such loading actions.
- Actifuse Microgranules should not be used in volumetrically unconstrained sites (so the graft material cannot move or escape).

Sterility:

- Sterilised by gamma irradiation.
- Do not resterilise.
- If the packaging or seal are damaged do not use. Discard.

Contraindications:

- Impaction grafting for failed total hip or knee arthroplasty
- Direct loading of graft material i.e. In the absence of conventional implants, such as screws and rods
- Infection
- Inability to cover or deliberate non-coverage of graft site using soft tissue
- Avascular or compromised vascular network sites
- Patient metabolism may compromise bony regeneration
- Avoid use in patients where in the surgeon's opinion patient lifestyle, compliance and/or physical attributes would compromise clinical outcome
- Medication which could slow bone healing

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