Actifuse®

Shape

GB	INSTRUCTIONS FOR USE.	p1
DE	GEBRAUCHSANLEITUNG.	p3
GR	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.	p5
ES	INSTRUCCIONES DE USO.	p7
NO	BRUKSANVISNING.	p9
PT	INSTRUÇÕES DE UTILIZAÇÃO.	p11
FI	KÄYTTÖOHJEET.	p13
IT	ISTRUZIONI PER L'USO.	p15
FR	MODE D'EMPLOI.	p17
DK	BRUGSVEJLEDNING.	p19
SE	BRUKSANVISNING.	p21
NL	GEbruIKSINSTRUCTIES.	p23
IL	הוראות לשימוש	p25
PO	INSTRUKCJA UŻYCIA	p27
CZ	POKYNY PRO POUŽITÍ.	p29
BU	УКАЗАНИЯ ЗА УПОТРЕБА.	p31
HU	HASZNÁLATI ÚTMUTATÓ.	p33
EST	KASUTUSJUHISED.	p35
LI	NAUDOJIMO INSTRUKCIJOS.	p37
RO	INSTRUCȚIUNI DE UTILIZARE.	p39
SK	NÁVOD NA POUŽITIE.	p41

Revision number 25012 R1

Apatech®
Redefining orthobiologics

INSTRUCTIONS FOR USE

1

Phase-pure, porous silicated calcium phosphate granules for use as a synthetic bone graft substitute.

Sterile - Single use - Store in a cool, dry place.

Description:

- 0.8 wt% Silicated calcium phosphate
- Total porosity 80% \pm 2.5%
- Water-soluble alkylene oxide copolymer carrier

Indications:

- Bone graft substitutes are intended to be used in place of corticocancellous, or cancellous allograft or autograft bone.
- The mechanical environment for such uses experience either low load requirements or compression.
- Typical surgical applications for bone graft substitutes are:-
 - Small void filling, e.g. after removal of a small bone tumour or following bone fracture reduction or in osteotomies and plastic surgery.
 - Spinal fusion, where a cage or screw fixation device is used to relieve the graft site from physiological loads.
- It is not intended to be used in place of cortical strut allograft bone where high tensile, torsion and/or bending strength are required. The products are used by orthopaedic surgeons in place of allograft bone (bone from humans stored in bone banks).

Directions for Use:

- Open both outer (non-sterile) and inner (sterile) pouches.
- Use in place of/or added to morsellised autograft bone.
- Use in place of/or added to morsellised allograft bone.
- Actifuse Shape is designed to be used alone. However, it may be mixed with sterile saline, autologous blood or bone marrow aspirate but this may have an effect on handling.
- Actifuse Shape implanted alone or when mixed is used in the same manner as morsellised allograft or autograft.

Risks and warnings:

- The risks with Actifuse Shape (e.g. infection, graft loosening, non-union, delayed union, repeat surgery) are the same as those for allograft bone. However, with Actifuse Shape there is no risk of disease transmission.

- Actifuse Shape should not be used where it could be subject to ² tension, torsion, compression, shear or bending. A conventional implant (e.g., screw, rod) can protect the graft from such loading actions.
- Actifuse Shape should not be used in volumetrically unconstrained sites (so the graft material cannot move or escape).

Sterility:

- Sterilised by Electron Beam irradiation.
- Do not resterilise.
- If the packaging or seal are damaged do not use. Discard.

Contraindications:

- Impaction grafting for failed total hip or knee arthroplasty
- Direct loading of graft material, i.e. in the absence of conventional implants, such as screws and rods
- Infection
- Inability to cover or deliberate non-coverage of graft site using soft tissue
- Avascular or compromised vascular network sites
- Patient metabolism may compromise bony regeneration
- Avoid use in patients where in the surgeon's opinion patient lifestyle, compliance and/or physical attributes would compromise clinical outcome
- Medication which could slow bone healing

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ApaTech Ltd., 370 Centennial Avenue,